


Instruction Manual for Zero Fog Goggles

I. Product description

- AURORA goggles can prevent liquid splash, saliva, water drop transmission, bacteria, flying debris, sand or dust. It can be used in hospitals, laboratories, mines, gardens, barns, construction sites and other places.
- Adopt electronic anti-fogging design, actively defog, remove water vapor and residual heat out of the goggles, without fogging or sultry heat, allowing you to enjoy super clear vision, helping special groups such as doctors and nurses and special places solve the problem of fogging with goggles that affects vision.
Important note: Please press the switch of the zero fog goggles before use to turn on the electronic defogging function.
- The electronic filter of AURORA goggles adopts HEPA H13 filtration, 99% effective filtration of smoke, bacteria, and air purification.
- The unique air protection design of AURORA goggles allows clean air to fill the eyes to form a safe air protection zone.
- The 180-degree ultra-wide view field design provides an all-round field of view, protects your eyes from splash particles from all sides.
- Long battery life, continuous power supply for 8 hours after fully charged.
- Designed for long time use, the filter element can be replaced.
- Ergonomically designed and equipped with a comfortable adjustable elastic headband, which can be adjusted according to personal needs. The soft material on the contact surface ensures a perfect fit to every head shape and is completely sealed between the face and the environment. Comfortable to wear, suitable for everyone.
- Made of crystal clear soft food-grade PVC, can be comfortably installed on most prescription glasses, comfortably fit the glasses and crystal clear material ensuring that you can see objects clearly.
- High scratch resistance and defogging functions. The special coating of the lens makes the goggles very scratch-resistant on the outside, and no fog inside. Polarized lenses protect your eyes from harmful ultraviolet rays.



Instruction Manual for Zero Fog Goggles


II. Parameters of Zero Fog Goggles

- Product name: Zero Fog Goggles
- Product model: AM1+
- Input voltage: DC5V
- Input power: 0.8W
- Battery parameters: 3.7V 600mAh, polymer battery
- Charging time: 2.5 hours (with DC5V 500mAh)
- Charger specifications: DC5V 500mA or above
- Working time: 8 hours
- Goggles size: 200×90×88mm
- Goggles weight: 105g
- Packing method: pouch bag + box
- Goggles color: light blue and transparent
- Lens: PC anti-fog lens
- Whether to include battery: Yes

III. Charging/Working and LED indicator

Sequence	Working mode	Indicator color	Comment
1	Not working	LED indicator not light up	
2	Charging	LED indicator in red color	
3	Fully charged	LED indicator in green color	
4	Working	LED indicator in green color	
5	Stop working	LED indicator not light up	

IV. Wearing method



V. Precautions

- Gently peel off the protective layer of the lens before the first time use.
 - Press the switch of goggles on, and the indicator light turns green, indicating that it starts to work.
 - Cleaning and disinfection: firstly use a wet towel or wet paper towel to remove dust, sweat and other impurities on the goggles, then disinfect the goggles with alcohol or a disinfectant wet paper towel, and dry as soon as possible after disinfection.
 - Be careful not to clean or disinfect it when the goggles are dry to prevent alcohol from damaging the coating of the lens.
 - Please avoid liquid from flowing into the USB port. If you need to rinse with water, please take out the electronic filter before rinse.
- Kind Reminders:
- If the product has been contaminated or the wearer leaves the severely infected area, it should be disposed specially and not reused.
 - The anti-fog coating of the lens will be affected and the anti-fog effect of the lens will be weakened after cleaned and repeated disinfected by alcohol and wipe the lens. This is a material characteristic and is a normal phenomenon.